


Safety First!

Using MagHealy the Right Way



Safety Precautions

- Keep away from children
- Avoid immersion in liquids
- Do not use when driving or using heavy machinery
- When using MagHealy with implanted or invasive medical devices, seek clarification from the manufacturers of those medical devices before use.



Maintenance & Warranty

- For defects, contact the manufacturer
- Any attempt to repair the device voids the warranty
- Equipment used with MagHealy must be in good condition
- No maintenance while equipment is in use

User Guidelines

- Intended for adults only
- Not for pregnant women
- Not for those with epilepsy or heart issues
- Not if you have a pacemaker
- Not for use in the head area
- Not for direct skin contact

Device Care

- Store & transport as per instructions
- Designed for indoor use only
- Avoid using in shielded areas (ICUs, airplanes)
- Only use supplied accessories and compatible phones
- When not in use, protect the device against damage by keeping it closed and away from:
Children | Dust & sunlight | Moisture | Pets & vermin
- Keep apart from other equipment.



Usage Tip
MagHealy may heat up when charging. Wait 10 minutes after charging before using the device.